



# Стандарты CONOSTAN для анализа масел и нефтепродуктов

**CONOSTAN**<sup>®</sup>  
Oil Analysis Standards

**CONOSTAN**<sup>®</sup>  
Oil Analysis Standards

IS A DIVISION OF

**SCP SCIENCE**



## CALIBRATION AND VERIFICATION STANDARDS

METALS IN OIL STANDARDS

SULFUR IN CRUDE OIL

FLASH POINT STANDARDS

TAN & TBN STANDARDS

*“The benchmark for accurate measurements”*

## Calibration & Verification Standards for the Analysis of:

- Lubricants: New and In-service
- Lubricant Additives
- Petroleum Products
- Organic Fluids / Materials



IS A DIVISION OF

**SCP SCIENCE**

[post@energolab.com](mailto:post@energolab.com)

## QUALITY COMMITMENT

OUR GOAL IS TO SERVE YOU BETTER THROUGH OUR CONTINUOUS IMPROVEMENT

SCP SCIENCE, the manufacturer of CONOSTAN Standards is very proud to announce that the company is registered/ accredited to three international Quality Management System standards: ISO 9001, ISO 17025 and ISO Guide 34.



This comprehensive coverage means that you will be receiving a quality product no matter what you purchase from SCP SCIENCE / CONOSTAN 's extensive offering, each and every time you order. We are committed to providing quality products and services upon which customers can depend.

Through our all-encompassing ISO 9001 certification and laboratory focused ISO 17025 / Guide 34 accreditations, SCP SCIENCE / CONOSTAN has proven that it can provide the best and most reliable products and reference materials available in the industry. At SCP SCIENCE / CONOSTAN , we are truly dedicated to providing and exceeding the quality you require. SCP SCIENCE is a member of the Chemical Reference Material Manufacturers Association (CRMMA); an association of the major manufacturers of calibration standards whose mandate is "...to promote the production and marketing of high quality materials and to promote continuous improvement of industry products and standards..."

SCP SCIENCE is a leader in the field of standards and volumetric solutions and yet one of our key goals is to continually improve to serve you better. Our websites include many new products in response to your requests and to the new requirements of governments and organizations. Should you require a product that is not listed on our websites, please do not hesitate to contact our Customer Service Group.



## TABLE OF CONTENTS

QUALITY COMMITMENT.....	3
ABOUT CONOSTAN .....	5
Our Products.....	5
METALLO-ORGANIC STANDARDS .....	6
Applications .....	6
Single-Element Standards .....	7
D-Series Standards for Joint Oil Analysis Program (JOAP) .....	9
Used Oil - Certified Reference Material - EnviroMAT™ .....	11
BIODIESEL STANDARDS.....	12
Metals in Biodiesel.....	12
Sulfur in Biodiesel .....	13
METALS IN OIL STANDARDS.....	14
Metals in Crude Oil.....	14
Metals in Residual Oil.....	14
PREMISOLV ICP SOLVENT.....	15
CHLORINE IN OIL STANDARDS.....	15
SULFUR IN OIL STANDARDS .....	16
Sulfur in Isooctane .....	16
Sulfur in Residual Oil .....	17
Sulfur in Crude Oil .....	17
Sulfur Standards in Mineral Oil and Diesel Fuel (100 g) .....	18
VISCOSITY STANDARDS .....	19
FLASH POINT STANDARDS .....	23
PARTISTAN - PARTICLE COUNTER STANDARDS.....	24
FTIR STANDARDS .....	24
TOTAL ACID NUMBER (TAN) STANDARDS .....	25
TOTAL BASE NUMBER (TAN) STANDARDS .....	26
CUSTOM QUOTATION REQUEST FORM FOR METALLO-ORGANIC STANDARDS .....	27

## ABOUT CONOSTAN

### OUR REPUTATION

SCP SCIENCE is proud to present CONOSTAN - the world's most trusted name in oil standards, whose industry-leading position is the result of a uniquely superior product chemistry and manufacturing technology together with proprietary blending techniques.

SCP SCIENCE, the manufacturer of CONOSTAN Standards is very proud to announce that the company is registered/ accredited to three international Quality Management System standards: ISO 9001, ISO 17025 and ISO Guide 34.

This comprehensive coverage means that you will be receiving a quality product no matter what you purchase from SCP SCIENCE / CONOSTAN 's extensive offering, each and every time you order. We are committed to providing quality products and services upon which customers can depend.

Through our all-encompassing ISO 9001 certification and laboratory focused ISO 17025 / Guide 34 accreditations, SCP SCIENCE/ CONOSTAN has proven that it can provide the best and most reliable products and reference materials available in the industry. At SCP SCIENCE / CONOSTAN , we are truly dedicated to providing and exceeding the quality you require.

SCP SCIENCE is a member of the Chemical Reference Material Manufacturers Association (CRMMA); an association of the major manufacturers of calibration standards whose mandate is "...to promote the production and marketing of high quality materials and to promote continuous improvement of industry products and standards..."

SCP SCIENCE is a leader in the field of standards and volumetric solutions and yet one of our key goals is to continually improve to serve you better. Our websites ([www.scpscience.com](http://www.scpscience.com) and [www.CONOSTAN.com](http://www.CONOSTAN.com)) include many new products in response to your requests and to the new requirements of governments and organizations. Should you require a product that is not listed on our websites, please do not hesitate to contact our Customer Service Group.



### OUR PRODUCTS

Oil standards are used extensively in the calibration and operation of instruments that analyze elements in oil and other organic fluids. Our product lines offer a wide variety of reference standards, solvents, and reagents for ICP, DCP, rotating disk electrode, XRF, AA, and other analytical spectrometric techniques.

#### Our products are optimized for:

- Compatibility – We offer combinations of 38 different elements over an extensive concentration range.
- Solubility – Our standards are soluble in a variety of substances: ketones, mineral oil, xylene, kerosene, etc. We also produce blank oils and PremiSolv, a kerosene alternative, for use as solvents.
- Volatility – Our standards are made in ultra-pure, highly processed hydrocarbon oil only: no solubilizers are used, making our standards extremely stable to volatile loss.
- Viscosity – The viscosity range of our standards at room temperature is ideal for instrumental applications.
- Instrumental response – Our standards provide excellent analytical response over a wide range of applications.
- Shelf life – All CONOSTAN element standards and spectroscopy products have a one-year minimum shelf life from the date of shipment.


## METALLO-ORGANIC STANDARDS

### APPLICATIONS

CONOSTAN standards are used extensively in a variety of industries: energy, environmental, aircraft, railroad, automotive, heavy equipment, mining, chemical and others. Essentially, wherever calibration of instruments for analyzing elements in oil and other organics is needed, CONOSTAN standards provide consistent composition and performance.

#### Typical Uses:

1. For ICP, DCP, AA and rotating disk emission spectrometric determinations of trace metals in organic materials, such as wear metals in used lubricating oils.
2. In AA and ICP spectrometric analyses of oils and other organics for As, Ba, Cd, Cr, Pb, Se, and other metals of environmental concern.
3. In AA and ICP emission analyses for trace metals in organics such as Co and Ni from metal-catalyzed reactions.
4. For XRF determinations of multi-element mixtures in organic systems such as Fe, Ni, Cu, and V in crude oils.
5. For emission spectrometric determinations of additive metals such as Ba, Ca, Mg, P, and Zn for quality control of motor oil formulations.
6. As internal standards and matrix adjustment components, and for the preparation of special multi-element blended standards to meet specific calibration requirements.



Eliminate the need to change ICP Sample Introduction components for different applications by using the

#### **MINI X-FLOW NEBULIZER**

Whether your samples are new or used oils, clean waters, suspensions, sludges or in HF, the Mini X-Flow Nebulizer (Cat. # 020-06.0-6.11) can handle them all - and with the best performance coupled to a Cyclonic Spray Chamber.

## SINGLE-ELEMENT STANDARDS

BLENDED IN 20 CST BLANK OIL (SIZE: 50 G)

Elements	1000 ppm	5000 ppm
	50 g	50 g
Ag	150-100-475	150-500-475
Al	150-100-135	150-500-135
As	150-100-335	-
B	150-100-055	150-500-055
Ba	150-100-565	150-500-565
Be	150-100-045	150-500-045
Bi	150-100-835	150-500-835
Ca	150-100-205	150-500-205
Ce	150-100-585	150-500-585
Cd	150-100-485	150-500-485
Co	150-100-275	150-500-275
Cr	150-100-245	150-500-245
Cu	150-100-295	150-500-295
Fe	150-100-265	150-500-265
Hg	150-100-805	-
In	150-100-495	150-500-495
K	150-100-195	150-500-195
La	150-100-575	150-500-575

Elements	1000 ppm	5000 ppm
	50 g	50 g
Li	150-100-035	150-500-035
Mg	150-100-125	150-500-125
Mn	150-100-255	150-500-255
Mo	150-100-425	150-500-425
Na	150-100-115	150-500-115
Ni	150-100-285	150-500-285
P	150-100-155	150-500-155
Pb	150-100-825	150-500-825
Sb	150-100-515	150-500-515
Si	150-100-145	150-500-145
Sn	150-100-505	150-500-505
Sr	150-100-385	150-500-385
Ti	150-100-225	150-500-225
V	150-100-235	150-500-235
W	150-100-745	150-500-745
Y	150-100-395	150-500-395
Zn	150-100-305	150-500-305

	100 ppm
Sc*	150-103-341
As*	150-103-331

\* size: 100 g, 20 cSt

Elements	1000 ppm	5000 ppm
As***	150-101-331	150-102-331
Hg***	150-101-801	150-101-805
Se***	150-101-341	-

\*\*\* size: 100 g, 75 cSt

	2000 ppm
Sc*	150-500-215

\*\* size: 50 g, 20 cSt



## MULTI-ELEMENT STANDARDS

BLENDED IN 75 CST BLANK OIL

See detailed Certificate of Analysis on page 29

ppm	S-12		S-21		S-21+K	
	100 g	200 g	100 g	200 g	100 g	200 g
10	150-012-001	150-012-009	150-021-002	150-021-018	150-021-042	150-021-051
30	150-012-004	150-012-012	150-021-008	150-021-027	150-021-045	150-021-056
50	150-012-006	150-012-014	150-021-010	150-021-030	150-021-047	150-021-058
100	150-012-002	150-012-010	150-021-003	150-021-019	150-021-043	150-021-052
300	150-012-005	150-012-013	150-021-009	150-021-028	150-021-046	150-021-037
500	150-012-007	150-012-015	150-021-011	150-021-031	150-021-048	150-021-059
900*	150-012-008	150-012-016	150-021-015	150-021-035	150-021-049	150-021-061

\* 900 ppm nominal value, 885 ppm actual value.

**S-12:** Ag, Al, Cr, Cu, Fe, Mg, Na, Ni, Pb, Si, Sn, Ti

**S-21:** Ag, Al, B, Ba, Ca, Cd, Cr, Cu, Fe, Mg, Mn, Mo, Na, Ni, P, Pb, Si, Sn, Ti, V, Zn

ppm	AM-Special		
	100 g	200 g	400 g
500	150-250-010	150-250-023	-
900	150-250-014	150-250-027	150-250-004
1000	150-250-006	150-250-016	150-250-001
2000	-	150-250-017	-
2500	150-250-007	150-250-018	-
3000	150-250-009	150-250-020	-
4000	-	150-250-021	-
5000	150-250-011	150-250-024	150-250-003
6000	150-250-012	-	-
7000	150-250-013	150-250-025	-

**AM-Special:** Ba, Ca, Mg, P, Zn

## CUSTOM BLENDS

If you require a combination or concentration of elements that we do not routinely stock, custom blends are available and can be shipped within two working days. To order, you may fill out and mail the form on page 27 of this booklet. We can make:

- Any combination of elements listed here at custom concentrations
- Single-element standards at custom concentrations
- S-21, S-12, and AM-Special blends at custom concentrations
- Additions to S-21, S-12, and AM-Special blends

Size: available in 100 g, 200 g, and 400 g



## METALLO-ORGANIC STANDARDS (CONT'D)

### D-SERIES STANDARDS FOR JOINT OIL ANALYSIS PROGRAM (JOAP)

Direct from the Original Source

CONOSTAN D-Series standards are available directly from SCP SCIENCE. Each D3, D12, and D19 standard is furnished with a Certificate of Analysis.

CONOSTAN is the original source of D-Series standards. More than 35 years ago, the U.S. Department of Defense's Spectrometric Oil Analysis Program Standards Committee required standards for its wear metals analysis programs.

With no reliable commercial source of metals in oil standards, CONOSTAN's research department set to work in developing a reliable standard. The result was CONOSTAN's uniquely superior product chemistry, which was adopted by the Department of Defense for its D-Series standards.

As with all CONOSTAN products, you are guaranteed that the D-Series of standards are extremely stable and accurate.



### CROSS REFERENCE WITH U.S. DEPARTMENT OF DEFENSE STOCK NUMBERS

CONOSTAN product	Dept. of Defense NSN	Catalogue #
<b>D3-100</b>	9150-01-283-0249	150-300-019
<b>D12-5</b>	9150-01-307-3343	150-300-005
<b>D12-10</b>	9150-00-179-5145	150-300-001
<b>D12-30</b>	9150-00-179-5144	150-300-003
<b>D12-50</b>	9150-00-179-5143	150-300-006
<b>D12-100</b>	9150-00-179-5142	150-300-002
<b>D12-300</b>	9150-00-179-5141	150-300-004
<b>D19-0</b>	9150-00-179-5137	150-300-008
<b>D19 set</b>	9150-01-355-1178	150-300-018

## METALLO-ORGANIC STANDARDS (CONT'D)

### D-SERIES STANDARDS

ppm	D3	D12	D19
	200 g	200 g	100 g
0 (Blank)	150-300-030	150-300-030	150-300-008
5	-	150-300-005	150-300-013
10	-	150-300-001	150-300-009
30	-	150-300-003	150-300-011
50	-	150-300-006	150-300-014
100	150-300-019	150-300-002	150-300-010
300	-	150-300-004	150-300-012
500	-	-	150-300-015
700	-	-	150-300-016
900	-	-	150-300-017

\* Custom ppm preparations available.

**D3:** B, Mo, Zn    **D12:** Ag, Al, Cr, Cu, Fe, Mg, Na, Ni, Pb, Si, Sn, Ti

**D19:** Ag, Al, B, Ba, Cd, Cr, Cu, Fe, Mg, Mn, Mo, Na, Ni, Pb, Si, Sn, Ti, V, Zn

**Blank oil for D-Series:** 150-300-008

### D19 SET

Product	Catalogue #
D19 set	150-300-018

The D19 set comprises the following concentrations and quantities:

ppm	Quantity	D19 (100 g)
0 (blank)	4	150-300-008
5	1	150-300-013
10	1	150-300-009
30	1	150-300-011
50	1	150-300-014
100	3	150-300-010
300	2	150-300-012
500	1	150-300-015
700	1	150-300-016
900	1	150-300-017



## METALLO-ORGANIC STANDARDS (CONT'D)

### CUSTOM BLENDS OF D-SERIES STANDARDS

If you require a combination or concentration of elements that we do not routinely stock, custom blends are available.

**We can make:**

- Any combination of elements listed previously as custom concentrations
- Single-element standards at custom concentrations
- D3, D12, and D19 blends at custom concentrations
- Additions to D3, D12, and D-19 blends

**Size:** Available in 100 g or 200 g

**Available elements:** Ag, Al, B, Ba, Be, Bi, Ca, Cd, Co, Cr, Cu, Fe, In, K, La, Li, Mg, Mn, Mo, Na, Ni, P, Pb, Sb, Si, Sn, Sr, Ti, V, W, Y, Zn

### USED OIL - CERTIFIED REFERENCE MATERIAL - ENVIROMAT™

EnviroMAT™ Certified Reference Materials (CRM) can be an invaluable component of any laboratory quality control program. Consensus certification removes any chance of analytical bias. Used as a method performance Quality Control standard and / or QC in oil analysis by ICP-OES and AA.

- Each CRM is certified through a round-robin study employing specific methods of analysis:
  - Independant verification from multiple laboratories
- Includes Certificate of Analysis listing consensus values, confidence and tolerance intervals, and instructions for use: Complete documentation for audit purposes

EnviroMAT™ Standars	Symbol	Quantity	Catalogue #
As***	HU-1	125 ml	140-025-041



## BIODIESEL STANDARDS

Biodiesel is a clean burning alternative fuel, produced from domestic, renewable resources, such as virgin soybean and rapeseed (canola) oils. Biodiesel contains no petroleum, but it can be blended at any concentration with petroleum diesel to create a biodiesel blend. Blends of biodiesel and pure unblended biodiesel can be used in modern diesel engines with no modification required.

SCP SCIENCE-CONOSTAN offers a new line of biodiesel standards for the analysis of metals and sulfur in biodiesel fuel. Manufactured in accordance to ASTM methods D7039, D6.751, D5453 and EN14214 for ICP and XRF analysis. Each Standard is supplied a Certificate of Analysis.



## METALS IN BIODIESEL

SCP SCIENCE-CONOSTAN offers multi-element standards for Ca, K, Mg, Na and P in B100 Biodiesel. Blank standard is also available. Manufactured in accordance to methods D4951, EN14531, EN14017, EN14108 and EN14109 for ICP analysis. Single element, multi element and custom multi element in biodiesel standards are available upon request.

METALS IN BIODIESEL	100% BIODIESEL (B100) BDM5 (100 g)	100% BIODIESEL (B100) BDM2A (100 g)	100% BIODIESEL (B100) BDM2B (100 g)
Elements	Ca, K, Mg, Na, P	K, Na	Ca, Mg
ppm	Catalogue #		
0 (Blank)	150-441-000	150-441-000	150-441-000
2.5	150-441-005	150-441-030	150-441-065
5	150-441-010	150-441-035	150-441-070
10	150-441-015	150-441-040	150-441-075
15	150-441-020	150-441-045	150-441-080
20	150-441-025	150-441-050	150-441-085
25	-	150-441-055	150-441-090
50	-	150-441-060	150-441-095



### SULFUR IN BIODIESEL

CONOSTAN offers sulfur in biodiesel standards in 5% (B5) and 20% (B20) biodiesel blends. Blank standard is also available. Manufactured in accordance to ASTM methods D7039, D6.751, D5453 and EN14214 for ICP and XRF analysis. Custom blended biodiesel standards are available upon request.

SULFUR IN BIODIESEL	5% BIODIESEL (B5) SULFUR (100 g)	20% BIODIESEL (B20) SULFUR (100 g)	100% BIODIESEL (B100) SULFUR (100 g)
ppm	Catalogue #		
0 (Blank)	150-440-000	150-440-050	150-440-100
5	150-440-005	150-440-055	150-440-105
10	150-440-010	150-440-060	150-440-110
15	150-440-015	150-440-065	150-440-115
30	150-440-020	150-440-070	150-440-120
50	150-440-025	150-440-075	150-440-125
75	150-440-030	150-440-080	150-440-130
100	150-440-035	150-440-085	150-440-135
200	150-440-040	150-440-090	150-440-140
500	150-440-045	150-440-095	150-440-145



## METALS IN OIL STANDARDS

### METALS IN CRUDE OIL

Metals in Crude Oil for use in accordance with ASTM D5708 and D586.3. Complete with a Certificate of Analysis.

Iron, ppm	Nickel, ppm	Complete Set	150-451-060
		Vanadium, ppm	Catalogue # (100 ml)
1	2	2	150-451-000
300	10	500	150-451-005
500	100	25	150-451-010
100	80	250	150-451-015
200	40	100	150-451-020
400	5	400	150-451-025
1	60	300	150-451-030
500	2	200	150-451-035
100	100	2	150-451-040
300	50	250	150-451-045
200	20	500	150-451-050
50	100	50	150-451-055

### METALS IN RESIDUAL OIL

Metals in Residual Oil for use in accordance with ASTM D5708 and D586.3. Complete with a Certificate of Analysis.

Iron, ppm	Nickel, ppm	Complete Set	150-421-060
		Vanadium, ppm	Catalogue # (100 ml)
1	2	2	150-421-000
300	10	500	150-421-005
500	100	25	150-421-010
100	80	250	150-421-015
200	40	100	150-421-020
400	5	400	150-421-025
1	60	300	150-421-030
500	2	200	150-421-035
100	100	2	150-421-040
300	50	250	150-421-045
200	20	500	150-421-050
50	100	50	150-421-055



## PREMISOLV ICP SOLVENT

PremiSolv is a zero-odor alternative to kerosene or xylene for use as a diluent or zero-point standard in ICP/DCP analysis of metals in oil and other organic fluids.

### PremiSolv Features:

- Extremely low odor – for a safer, more comfortable working environment
- Extremely low toxicity – compared with kerosene or xylene
- Extremely low metal content – comes with a certificate of analysis listing the concentrations of 34 different metals including sulfur
- Non-hazardous – for shipping and disposal



Size	Catalogue #
3.78 L (1 gallon)	150-700-003
18.9 L (5 gallons)	150-700-002

**WANT TO SEE THE DIFFERENCE YOURSELF? Ask for a free 400 ml sample. Cat. No. 150-700-000**

## CHLORINE IN OIL STANDARDS

The chlorine product line is designed for calibration of XRF, ICP, and other analytical techniques. Stocked concentrations are shown in the table below.

Custom concentrations are available.

### Chlorine Standards (size: 100 g)

ppm	Catalogue #
o (Blank)	150-200-008
10	150-200-001
100	150-200-002
500	150-200-005
1000	150-200-003
5000	150-200-006
10000	150-200-004
50000	150-200-007



## SULFUR IN OIL STANDARDS

CONOSTAN offers an expanding line of Sulfur standards that today encompasses several matrices described below and on the following pages.

The quality of many petroleum products are related to the amount of sulfur present. In fuels, the amount of sulfur is related to its performance characteristics and vehicle emission levels. Levels of sulfur in petroleum and petroleum products are regulated through federal, state and local agencies.



### SULFUR IN ISOCTANE

Sulfur in Isooctane by UV Fluorescence

ppm	Description	Catalogue #
<b>0, 0.5, 1.0, 2.5, 5.0, 7.5, 10</b>	For very low Sulfur level, set of 7 x 10 ml vials. Complete with a Certificate of Analysis.	150-430-010
<b>0, 5.0, 10, 25, 50, 100, 250</b>	For low Sulfur level, set of 7 x 10 ml vials. Complete with a Certificate of Analysis.	150-430-020
<b>0, 50, 100, 250, 500, 750, 1000</b>	For high Sulfur level, set of 7 x 10 ml vials. Complete with a Certificate of Analysis.	150-430-030

### SULFUR IN ISOCTANE

Sulfur in Isooctane by X-Ray Fluorescence

ppm	Catalogue #
<b>o (Blank)</b>	150-430-101
<b>50</b>	150-430-102
<b>100</b>	150-430-103
<b>250</b>	150-430-104
<b>500</b>	150-430-105
<b>750</b>	150-430-106
<b>1000</b>	150-430-107
<b>Set of 7</b>	150-430-100





## SULFUR IN OIL STANDARDS (CONT'D)

### SULFUR IN RESIDUAL OIL

ppm	Catalogue # (50 ml)	Catalogue # (100 ml)
2500	150-420-100	150-420-005
3500	150-420-105	150-420-010
5000	150-420-110	150-420-015
7500	150-420-120	150-420-020
10000	150-420-125	150-420-025
15000	150-420-130	150-420-030
20000	150-420-135	150-420-035
25000	150-420-140	150-420-040
30000	150-420-145	150-420-045
35000	150-420-150	150-420-050
40000	150-420-155	150-420-055
50000	150-420-160	150-420-060

### SULFUR IN CRUDE OIL

We have 10 standards from 500 to 50000 ppm. Sulfur in Crude Oil is in accordance with Method ASTM D26.22 and D4294. Complete with a Certificate of Analysis.

ppm	Catalogue # (100 ml)
500	150-450-000
1,000	150-450-005
1,500	150-450-020
2,500	150-450-010
5,000	150-450-015
10,000	150-450-025
20,000	150-450-030
30,000	150-450-035
40,000	150-450-040
50,000	150-450-045





## SULFUR IN OIL STANDARDS (CONT'D)

### SULFUR STANDARDS IN MINERAL OIL AND DIESEL FUEL (100 G)

#### Sulfur in Mineral Oil

The sulfur in oil product line is designed for calibration of XRF, ICP, and other analytical instruments according to various ASTM methods (such as D26.22, D3246., D4294, D5453, D6.334, and D6.443). Stocked concentrations are shown in the table below.

#### Sulfur in Diesel Fuel

Our line of sulfur in diesel fuel is specially engineered to have an elevated flash point, making it suit-able for shipping as a non-hazardous product. Stocked concentrations are shown in the table below.

#### Custom Concentrations

For concentrations not listed in the table below, custom concentrations are available in both mineral oil and diesel fuel. Contact us for further information.

ppm	In Mineral Oil (100 g)	In Diesel Fuel (100 g)
0 (Blank)	150-400-025	150-410-012
5	CB2-000-018	150-410-008
10	150-400-001	150-410-001
25	150-400-009	150-410-013
50	150-400-018	150-410-009
100	150-400-002	150-410-002
250	150-400-010	CB2-000-023
500	150-400-019	150-410-010
750	150-400-023	150-410-018
1000	150-400-003	150-410-003
2500	150-400-011	CB2-000-024
5000	150-400-020	150-410-011
7500	150-400-024	150-410-021
10000	150-400-004	150-410-004
15000	150-400-005	150-410-006
20000	150-400-008	150-410-007
25000	150-400-012	CB2-000-028
30000	150-400-014	150-410-022
40000	150-400-016	150-410-023
50000	150-400-021	150-410-024

\* Custom ppm preparations available

\*\* A Certificate of Analysis reporting concentration and density are included with each standard

## VISCOSITY STANDARDS

The same product reliability and stability that our customers trust in metallo-organic standards are now available in General Purpose Viscosity standards. These certified, mineral oil based, viscosity standards were developed for calibration and verification of all types of viscometers, such as glass capillary viscometers.

All standards are traceable to National Standards in accordance to ASTM & IP methods. Each standard carries a two year stability guarantee. The determination of kinematic and dynamic viscosity were made in accordance with ASTM D445/446. and ISO 3104/3105, ISO/IEC 17025 and are traceable to the NIST (National Institute of Standards and Technology).



See detailed Certificate of Analysis on page 28

Viscosity Standards	VOLUME				
	125 ml	500 ml	1 Litre	4 Litre	20 Litre
<b>S3</b>	150-600-351	150-600-352	150-600-353	150-600-354	150-600-355
<b>S6</b>	150-600-141	150-600-142	150-600-143	150-600-144	150-600-145
<b>N4</b>	150-600-441	150-600-442	150-600-443	150-600-444	150-600-445
<b>N10</b>	150-600-181	150-600-182	150-600-183	150-600-184	150-600-185
<b>S20</b>	150-600-221	150-600-222	150-600-223	150-600-224	150-600-225
<b>N35</b>	150-600-261	150-600-262	150-600-263	150-600-264	150-600-265
<b>N44</b>	150-600-461	150-600-462	150-600-463	150-600-464	150-600-465
<b>S60</b>	150-600-301	150-600-302	150-600-303	150-600-304	150-600-305
<b>N100</b>	150-600-341	150-600-342	150-600-343	150-600-344	150-600-345
<b>S200</b>	150-600-231	150-600-232	150-600-233	150-600-234	150-600-235
<b>N350</b>	150-600-361	150-600-362	150-600-363	150-600-364	150-600-365
<b>N415</b>	150-600-471	150-600-472	150-600-473	150-600-474	150-600-475
<b>S600</b>	150-600-241	150-600-242	150-600-243	150-600-244	150-600-245
<b>N1000</b>	150-600-371	150-600-372	150-600-373	150-600-374	150-600-375
<b>S2000</b>	150-600-381	150-600-382	150-600-383	150-600-384	150-600-385
<b>N4000</b>	150-600-391	150-600-392	150-600-393	150-600-394	150-600-395
<b>S8000</b>	150-600-401	150-600-402	150-600-403	150-600-404	150-600-405
<b>N15000</b>	150-600-411	150-600-412	150-600-413	150-600-414	150-600-415
<b>S30000</b>	150-600-421	150-600-422	150-600-423	150-600-424	150-600-425



## VISCOSITY STANDARDS (CONT'D)

	KINEMATIC VISCOSITY IN MM <sup>2</sup> /S (CENTISTOKES)									Saybolt Viscosity
	20 °C/ 68 °F	25 °C/ 77 °F	37.78 °C/ 100 °F	40 °C/ 104 °F	50 °C/ 122 °F	60 °C/ 140 °F	80 °C/ 176 °F	98.89 °C/ 210 °F	100 °C/ 212 °F	37 °C/ 100 °F
<b>S3</b>	4.5	4.0	3.0	2.8	2.4	2.0	1.5	1.2	1.2	
<b>S6</b>	10	8.8	6.0	5.7	4.5	3.6	2.5	1.9	1.9	
<b>N4</b>	6.7	5.8	4.2	4.0	3.2	2.6	1.9	1.5	1.4	
<b>N10</b>	21	17	11	10	7.5	5.8	3.7	2.7	2.6	
<b>S20</b>	46	35	20	18	13	9.0	5.6	3.6	3.5	87
<b>N35</b>	90	67	36	32	21	15	8.4	5.4	5.3	167
<b>N44</b>	110	86	48	44	30	21	12	7.6	7.4	220
<b>S60</b>	160	119	60	54	35	26	12	7.7	7.5	281
<b>N100</b>	318	228	110	97	60	39	20	11	11	509
<b>S200</b>	715	487	206	180	103	64	30	17	16	954
<b>N350</b>	1400	940	370	330	180	110	46	24	23	1730
<b>N415</b>	1900	1200	480	410	220	130	55	29	28	2200
<b>S600</b>	2400	1600	600	520	280	160	66	34	32	
<b>N1000</b>	5100	3300	1200	1000	520	290	110	52	50	
<b>S2000</b>	8200	5200	1900	1600	780	400	150	70	68	
<b>N4000</b>	18000	11000	3900	3300	1600	840	280	123	117	
<b>S8000</b>	37000	23000	7900	6700	3200	1600	520	210	200	
<b>N15000</b>	64000	40000	13000	11000	5300	2700	850	340	320	
<b>S30000</b>		80000	28000	23000	11000	5800	1700	670	640	

All calibrations and tests are based on a master viscometer procedure located in ASTM D216.2 and the National Institute of Standards and Technology (NIST) value of 1.0016. mPa.s (centipoise) for water at 20°C (6.8°F). Custom standards are available. Please contact us for more information.



## VISCOSITY STANDARDS (CONT'D)

	DYNAMIC VISCOSITY IN MPA.S (CENTIPOISE)									Saybolt Viscosity
	20 °C/ 68 °F	25 °C/ 77 °F	37.78°C/ 100 °F	40 °C/ 104 °F	50 °C/ 122 °F	60 °C/ 140 °F	80 °C/ 176 °F	98.89°C/ 210 °F	100 °C/ 212 °F	37 °C/ 100 °F
<b>S3</b>	3.7	3.3	2.4	2.3	1.9	1.6	1.2	0.9	0.9	
<b>S6</b>	8.7	7.3	5.0	4.7	3.6	2.9	2.0	1.5	1.4	
<b>N4</b>	5.6	4.8	3.4	3.2	2.6	2.1	1.5	1.1	1.1	
<b>N10</b>	18	14	9.0	8.4	6.2	4.7	3.0	2.1	2.1	
<b>S20</b>	40	30	17	15	11	7.6	4.7	2.9	2.9	87
<b>N35</b>	78	59	31	28	18	13	7.0	4.4	4.3	167
<b>N44</b>	91	71	39	36	24	17	9.4	6.0	5.8	220
<b>S60</b>	138	102	52	46	30	22	9.9	6.3	6.1	281
<b>N100</b>	276	197	94	83	51	33	16	9.4	9.4	509
<b>S200</b>	613	416	174	152	87	54	24	15	13	954
<b>N350</b>	1200	810	320	280	150	92	38	20	19	1730
<b>N415</b>	1600	1100	410	350	190	110	45	23	23	2200
<b>S600</b>	2100	1400	510	440	240	140	55	28	26	
<b>N1000</b>	4400	2400	1000	940	440	240	92	43	41	
<b>S2000</b>	7200	4500	1600	1400	670	340	130	58	56	
<b>N4000</b>	16000	9700	3400	2900	1400	720	240	100	98	
<b>S8000</b>	33000	20000	6900	5900	2800	1400	440	180	170	
<b>N15000</b>	57000	36000	11000	9700	4700	2400	730	290	270	
<b>S30000</b>		72000	25000	20000	9700	5100	1500	570	550	



## VISCOSITY STANDARDS (CONT'D)

	DENSITY IN G/ML IN ACCORDANCE WITH ASTM D7042								
	20 °C/ 68 °F	25 °C/ 77 °F	37.78°C/ 100 °F	40 °C/ 104 °F	50 °C/ 122 °F	60 °C/ 140 °F	80 °C/ 176 °F	98.89°C/ 210 °F	100 °C/ 212 °F
<b>S3</b>	0.833	0.830	0.821	0.820	0.813	0.806	0.792	0.779	0.778
<b>S6</b>	0.831	0.828	0.819	0.818	0.811	0.805	0.791	0.779	0.778
<b>N4</b>	0.834	0.831	0.822	0.820	0.814	0.807	0.793	0.780	0.780
<b>N10</b>	0.842	0.839	0.831	0.829	0.823	0.816	0.804	0.791	0.791
<b>S20</b>	0.871	0.868	0.860	0.859	0.852	0.846	0.833	0.821	0.820
<b>N35</b>	0.872	0.869	0.861	0.860	0.853	0.847	0.834	0.823	0.822
<b>N44</b>	0.828	0.825	0.817	0.816	0.809	0.803	0.791	0.779	0.779
<b>S60</b>	0.863	0.860	0.852	0.851	0.845	0.839	0.826	0.815	0.814
<b>N100</b>	0.867	0.864	0.857	0.855	0.849	0.843	0.831	0.820	0.819
<b>S200</b>	0.858	0.855	0.847	0.846	0.840	0.834	0.822	0.910	0.810
<b>N350</b>	0.863	0.860	0.852	0.851	0.845	0.839	0.827	0.816	0.815
<b>N415</b>	0.865	0.862	0.854	0.853	0.847	0.841	0.830	0.819	0.818
<b>S600</b>	0.866	0.864	0.856	0.855	0.849	0.843	0.832	0.820	0.820
<b>N1000</b>	0.872	0.869	0.862	0.860	0.855	0.849	0.837	0.826	0.826
<b>S2000</b>	0.876	0.873	0.865	0.864	0.858	0.853	0.841	0.831	0.830
<b>N4000</b>	0.882	0.879	0.872	0.871	0.865	0.860	0.849	0.838	0.838
<b>S8000</b>	0.888	0.885	0.878	0.877	0.872	0.866	0.855	0.845	0.845
<b>N15000</b>	0.893	0.891	0.884	0.883	0.877	0.872	0.861	0.851	0.851
<b>S30000</b>		0.896	0.889	0.888	0.883	0.877	0.867	0.857	0.857



## FLASH POINT STANDARDS

Flash Point Standards for use in accordance with ASTM D93. Complete with a detailed Certificate of Analysis. Contact your local Sales Representative for volume discounts.

ASTM D 93 Nominal Flash Point (°C)	Catalogue # (80 ml)	Catalogue # (250 ml)
53	150-900-000	150-900-050
69	150-900-005	150-900-055
81	150-900-010	150-900-060
104	150-900-015	150-900-065
112	150-900-020	150-900-070
116	150-900-025	150-900-075
134	150-900-030	150-900-080
186	150-900-035	150-900-085
231	150-900-040	150-900-090
260	150-900-110	150-900-115



### TORCHES WITH FIXED ALUMINA INJECTOR

A 99.6% alumina injector is fused into a torch body offering a longer torch life. These models are compatible with Agilent®, Spectro® EOP and PerkinElmer®. Patent pending.





## PARTISTAN - PARTICLE COUNTER STANDARDS

PartiStan particle standards are designed for calibration and verification of automatic particle counters.

### Product History

In 1999, a new calibration procedure (ISO 11171) for automatic particle counters was introduced, rendering previous procedures (i.e. ISO 4402) obsolete. With the new procedure, primary calibration requires NIST SRM 2806. - a suspension of 3.3 mg/L of ISO medium test dust in super-clean hydraulic fluid. CONOSTAN offers a secondary standard for the cost-effective calibration and verification of automatic particle counters. PartiStan secondary standards are compliant with ISO 11171 and directly traceable to NIST SRM 2806.

Description	Size	Catalogue #
PartiStan 2806	400 ml	150-701-001
PartiStan resolution standard	400 ml	150-701-002
PartiStan SCF (super-clean fluid)	400 ml	150-701-003
PartiStan SCF (super-clean fluid)	3.78 L (1 gallon)	150-701-004
PartiStan UFTD (ultra-fine test-dust suspension in SCF)	400 ml	150-701-005

## FTIR STANDARDS

Our FTIR operational test standard is a petroleum oil-based fluid that looks and handles like routinely tested used-oil samples. It is designed for validating FTIR instrument performance in order to ensure repeatability and reproducibility.

Description	Size	Catalogue #
FTIR operational test standard	100 g	150-702-001


For initial setup, we include a disk containing the analytical methods necessary to perform tests with information specific to your instrument:

- DigiLab (Varian®) – available
- PerkinElmer® – not yet available (contact the manufacturer)
- Thermo Nicolet® – not yet available (contact the manufacturer)



## TOTAL ACID NUMBER (TAN) STANDARDS

The titration curves for every concentration of the CONOSTAN family of TAN standards have more prominent inflection points and a wider difference between the starting pH and final pH during the titration, in comparison with other suppliers. At minimum concentration for TAN standards, the pH change in the titration curve is now greater than even the highest concentrations of competitive standards. Because the inflection point is more prominent, CONOSTAN TAN standards offer uncertainties of less than 1% of stated concentration providing customers the opportunity for better calibration of their titrator and TAN titration method.



Description	Size	Catalogue #
0.05 mg/g KOH *	75 g	150-800-005
0.1 mg/g KOH *	75 g	150-800-011
0.5 mg/g KOH *	75 g	150-800-051
1.0 mg/g KOH *	75 g	150-800-101
1.5 mg/g KOH	75 g	150-800-151
2.0 mg/g KOH	75 g	150-800-205
2.5 mg/g KOH	75 g	150-800-255
3.0 mg/g KOH	75 g	150-800-305
4.5 mg/g KOH	75 g	150-800-455

\* Available in packs of 3.

## TOTAL BASE NUMBER (TAN) STANDARDS

Each CONOSTAN TBN standard is formulated to provide a strong leveraging influence on the titration curve creating dramatic and easy to determine inflection points. Concentration uncertainties for TBN are less than 1%.

Description	Size	Catalogue #
1.0 mg/g KOH	75 g	150-801-011
3.0 mg/g KOH	75 g	150-801-031
6.0 mg/g KOH	75 g	150-801-065
10.0 mg/g KOH	75 g	150-801-105
15.0 mg/g KOH	75 g	150-801-155
30.0 mg/g KOH	75 g	150-801-305
40.0 mg/g KOH	75 g	150-801-405
70.0 mg/g KOH	75 g	150-801-705



**CUSTOM MULTI-ELEMENT STANDARDS**  
**ORDER ONLINE**

SCP SCIENCE / CONOSTAN ISO 9001 Manufacturing Division can prepare custom multi-element calibration standards to exact specifications. Indicate the desired elements and their concentrations, together with the desired matrix, and we will recommend a mix stability with a minimum of 12 months. A Certificate of Analysis, reporting the actual values and not simply the calculated values, accompanies each standard.



## CUSTOM QUOTATION REQUEST FORM FOR METALLO-ORGANIC STANDARDS

Complete this form to receive a quotation for your specific oil-based standard. Photocopy for use with multiple requests.

### Contact Information

Name: \_\_\_\_\_

Company: \_\_\_\_\_

Mailing address: \_\_\_\_\_

City: \_\_\_\_\_ State/Province: \_\_\_\_\_

Country: \_\_\_\_\_ Zip/Postal code: \_\_\_\_\_

Telephone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_ Account number: \_\_\_\_\_

Please indicate the concentration (ppm) required for each element:

Elements	ppm	Elements	ppm	Elements	ppm
Ag		Fe		Pb	
Al		Hg		Sb	
As		In		Sc	
B		K		Se	
Ba		La		Si	
Be		Li		Sn	
Bi		Mg		Sr	
Ca		Mn		Ti	
Cd		Mo		V	
Co		Na		W	
Cr		Ni		Y	
Cu		P		Zn	

Size (g): \_\_\_\_\_

Special requirements: \_\_\_\_\_

Application: \_\_\_\_\_


Base Oil:  20cSt  75 cSt  245 cSt

Fax this form to: [post@energolab.com](mailto:post@energolab.com)



**CONOSTAN**<sup>®</sup>  
Oil Analysis Standards

Manufactured by SCP SCIENCE  
21800 Clark Graham, Suite 6110  
Quebec, Canada H9C 4S6  
Tel: 1-514-657-0261 / 1-(800) 361-6826



## Certificate of Analysis

**1.0 DESCRIPTION: CONOSTAN Viscosity Reference Standard, N350**  
 Catalogue Number: 150-600-361 / 150-600-362 / 150-600-363/ 150-600-364 / 150-600-365  
 Lot Number: 1015  
 Matrix: White mineral oil  
 Expiration Date: September 2015

**2.0 CERTIFIED VALUES AND ASSOCIATED UNCERTAINTY:**

Temperature		Kinematic Viscosity mm <sup>2</sup> /s (cSt)	Dynamic Viscosity mPa-s (cP)	Density (g/ml)	Saybolt Viscosity (SUS)
°C	°F				
20.00	68.00	1345	1160	0.8622	
25.00	77.00	888.8	763.5	0.8591	
37.78	100.00	350.5	298.5	0.8516	1624
40.00	104.00	303.1	257.7	0.8502	
50.00	122.00	166.3	140.4	0.8444	
60.00	140.00	99.09	83.09	0.8385	
80.00	176.00	42.27	34.94	0.8266	
98.89	210.00	22.54	18.38	0.8153	109.4
100.00	212.00	21.82	17.78	0.8147	

**\*Expanded Uncertainty**

Viscosity Range	Kinematic Viscosity mm <sup>2</sup> /s (cSt)	Dynamic Viscosity mPa-s (cP)
0.3 to 10	±0.12%	±0.12%
10 to 100	±0.17%	±0.17%
100 to 1000	±0.21%	±0.21%
1000 to 6000	±0.22%	±0.22%

Method of analysis and traceability:

This viscosity standard has been prepared according to ASTM methods D445, D446 and corresponding ISO methods 3104 and 3105. Kinematic viscosities have been determined using Master Viscometers calibrated according to ASTM method D2162 and based on the established kinematic viscosity of 1.0034mm<sup>2</sup>/s for distilled water at 20.00°C per ISO/TR3666. Conversion of Kinematic to Saybolt viscosity has been calculated according to ASTM method D2161. Density has been determined according to ASTM method D7042. Thermometers used for temperature measurements are NIST traceable.

\* The uncertainty of the certified values have been calculated from applicable uncertainty contributors (u). The combined uncertainty ( $u_c = \sqrt{\sum u_i^2}$ ) has been multiplied by a coverage factor (k) of 2 to provide a 95% confidence interval.

**3.0 REFERENCE VALUES:** None

**4.0 APPROVAL AND DATE OF CERTIFICATION:**

Certification Approval:

Certification Date:

Aiketa Mishra, Conostan Production Manager  
September 16, 2013



**CONOSTAN**<sup>®</sup>  
Oil Analysis Standards

Manufactured by SCP SCIENCE  
21800 Clark Graham, Bois d'Amie  
Quebec, Canada H9X 4B6  
Tel: 1-514-457-0701 / 1-(800) 361-4820

**Certificate of Analysis**



**1.0 DESCRIPTION:** CONOSTAN Multi-Element Standard  
S-21:10 ppm  
Catalogue Number: 150-021-002 / 150-021-018  
Lot Number: 21426010  
Matrix: Base Oil 75 cSt  
Expiration Date: 12 months from date of shipment (see bottle for date of shipment)

**2.0 CERTIFIED VALUES AND ASSOCIATED UNCERTAINTY:**

Certified Concentrations, ppm (µg/g):

Ag	10.00±0.04	Al	10.00±0.04	B	10.00±0.02	Ba	10.00±0.01
Ca	10.00±0.04	Cd	10.00±0.05	Cr	10.00±0.05	Cu	10.00±0.02
Fe	10.00±0.05	Mg	10.00±0.04	Mn	10.00±0.02	Mo	10.00±0.01
Na	10.00±0.03	Ni	10.00±0.04	P	10.00±0.05	Pb	10.00±0.08
Si	10.00±0.02	Sn	10.00±0.06	Ti	10.00±0.03	V	10.00±0.04
Zn	10.00±0.06						

Method of analysis and traceability:

This standard was prepared by weight measurements originating from assayed element Concentrates. A precursor blend was verified by atomic emission or absorption spectroscopy. Element concentrations for this standard are based on the Concentrate assay\* values and were prepared to within the uncertainty values listed above at the 95% Confidence Interval, as determined by weight measurements of blend components conducted on balances calibrated and verified with NIST traceable weights.

\*Each element Concentrate was assayed by classical wet chemical methods. Precision of assay measurement is ±0.5 percent maximum, but typically ±0.3 percent, or less. Assay accuracy is within one percent of measured value, but typically much less, as determined by co-measurement of, and traceability to, NIST Standards, or Certified Analytical Reagent Grade Chemicals, if no suitable NIST standards exists.

**3.0 REFERENCE VALUES:**

None

**4.0 APPROVAL AND DATE OF CERTIFICATION:**

Certification Date: January 15, 2014

Certification Approval:

Alketa Mbha, Chemist  
Conostan Production Manager



The CONOSTAN Division of SCP SCIENCE is the world's leading manufacturer of Elements in Oil standards and related spectroscopy products. CONOSTAN`s industry leading position is derived from superior product chemistry, manufacturing technology, and blending techniques. CONOSTAN products are used extensively in the calibration and operation of analytical instruments for the analysis of Elements in Oil and other organic fluids.



ВЫБИРАЯ НАС -  
ВЫ ВЫБИРАЕТЕ ЛУЧШЕЕ!

ООО «ЭНЕРГОЛАБ»

### ЦЕНТРАЛЬНЫЙ ОФИС

127018, г. Москва, 3-й проезд Марьиной Рощи,  
д. 40, к. 6, стр. 1

Тел: +7 (495) 926-02-15

Офис продаж:

e-mail: [post@energolab.com](mailto:post@energolab.com)

Сервисная служба:

[service@energolab.com](mailto:service@energolab.com)

### ПРЕДСТАВИТЕЛЬСТВО

660021, г. Красноярск,  
ул. Дубровинского, д. 110

Тел: +7 (391) 274-65-56



**CONOSTAN<sup>®</sup>**  
*Oil Analysis Standards*

[WWW.ENERGOLAB-AE.COM](http://WWW.ENERGOLAB-AE.COM)

